6065 TumoTrol, Level 1  
6066 TumoTrol, Level 2

4 x 2.0 ml  
4 x 2.0 ml

IVD For In Vitro Diagnostic Use  
In Vitro Diagnostika

## Intended Use

The TumoTrol, Level 1/2 are used as assayed control sera for Internal Quality Control (IQC) in the medical laboratory.<sup>1)2)</sup>

## Principle

The use of control materials is an objective assessment to monitor precision and accuracy and is an integral part of Good Laboratory Practices (GLP). Two levels of control are available to allow statistical quality control within the clinical range.

## Reagent

The TumoTrol, Level 1/2 are lyophilized control materials based on human plasma, defibrinated for the following analytes:

AFP, Beta-2-Microglobulin, CA 15-3, CA 19-9, CA 125, CA 72-4, CEA, Cyfra 21-1, Ferritin, hCG, PSA, f-PSA.

## Target Values and Ranges

The Target values are derived from replicate analysis in several laboratories and are specific for the lot of TumoTrol, Level 1/2. Ranges are calculated as target value  $\pm$  the maximal allowed deviation of the single measured value.<sup>2)</sup>

Target values and ranges are only valid for the indicated method and may slightly differ due to various test kits, calibrators and system applications.

## Procedure

The reconstituted TumoTrol, Level 1/2 should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument or reagent kit being used.

## Reconstitution

Remove screw cap from the vial and pipette exactly 2,0 ml dest. water to the lyophilisate.

Keep the vial for about 30 min at rest in a light-protected place until the lyophilisate has completely dissolved.

Swing carefully- do not shake! Avoid developing of foam.

## Storage and Stability

The TumoTrol, Level 1/2 are stable until the expiration date when stored unopened at 2-8 °C

After reconstitution: **8 hours at 18 – 25 °C**

**5 days at 2 – 8 °C**

**1 month at – 20 °C**

## Warnings and Precautions

The blood donations used for production are tested with CE-marked test kits and found to be non-reactive for HBsAg, anti-HIV-1/2 and anti-HCV. In addition HIV and HCV is tested by PCR. Despite of that the danger of infection can not be excluded with certainty. The control materials have to be treated with the same care as patient samples.<sup>3)</sup>

## Waste Management

Dispose of any discarded materials in accordance with the requirements of local regulations authorities.

## Limitations

The products should not be used past the expiration date.

The products are not intended for use as a calibrator.

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DiaServe Laboratories GmbH
Scil Diagnostics GmbH