6071 ElphoTrol, Level 1
6072 ElphoTrol, Level 2

4 x 2.0 ml
4 x 2.0 ml

For In Vitro Diagnostic Use
In Vitro Diagnostika

Intended Use

ElphoTrol, Level 1/2 are used as assayed control sera for Internal Quality Control (IQC) in medical laboratories¹⁾.

Principle

The use of control materials is an objective assessment to monitor precision and accuracy. Two levels of control are available to allow performance.

Reagent

ElphoTrol, Level 1/2 are liquid control materials based on human serum for the following analytes:

Albumin, α_1 -Globulin, α_2 -Globulin, β -Globulin, γ -Globulin.

Assigned Values and Control Ranges

Values and control ranges are based on replicate measurements in several laboratories and are specific for the lots of ElphoTrol, Level 1/2.

Assigned values are only valid for the indicated method and may differ due to various system applications.

Procedure

ElphoTrol, Level 1/2 should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument being used.

Reconstitution

After opening the vial ElphoTrol, Level 1/2 is ready for use.

Storage and Stability

ElphoTrol, Level 1/2 are stable until the expiration date when stored unopened at +2°/+8 °C.

After opening: **8 hours at +18°/+25 °C**
1 week at +2°/+8 °C

Warnings and Precautions

The blood donations used for production were tested with CE-marked test kits and found to be non-reactive for HBsAg, anti-HIV-1/2 und anti-HCV. In addition HIV and HCV were tested by PCR. Despite of that the danger of infection from biological material can not be excluded with certainty. The control materials have to be treated with the same care as patient samples²⁾.

Waste Management

Dispose of any discarded materials in accordance with the requirements of local regulations authorities.

Limitations

The products should not be used past the expiration date.

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DiaServe Laboratories GmbH
Scil Diagnostics GmbH