Hemoglobin A1c Controls

Hämoglobin A1c Kontrollen

REF

6079 HbA1cTrol, Level 1
6080 HbA1cTrol, Level 2

4 x 1.0 ml
4 x 1.0 ml

IVD For In Vitro Diagnostic Use
In Vitro Diagnostika

Intended Use

The HbA1cTrol, Level 1/2 are used as assayed control sera for Internal Quality Control (IQC) in the medical laboratory.¹⁾²⁾

Principle

The use of control materials is an objective assessment to monitor precision and accuracy and is an integral part of Good Laboratory Practices (GLP). Two levels of control are available to allow statistical quality control within the clinical range.

Reagent

The HbA1cTrol, Level 1/2 are lyophilized control materials based on human blood.

Target Values and Ranges

Target values are derived from replicate analysis in several laboratories and are specific for the lot of HbA1cTrol, Level 1/2. Ranges are calculated as target value \pm the maximal allowed deviation of the single measured value.²⁾

Target values and ranges are only valid for the indicated method and may slightly differ due to various test kits, calibrators and system applications.

Procedure

The reconstituted HbA1cTrol, Level 1/2 should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument or reagent kit being used.

Reconstitution

Remove screw cap from the vial and pipette exactly 1,0 ml dest. water to the lyophilisate.

Keep the vial for about 30 min at rest in a light-protected place until the lyophilisate has completely dissolved.

Swing carefully- do not shake! Avoid developing of foam.

Storage and Stability

The Hb A1cTrol, Level 1/2 are stable until the expiration date when stored unopened at 2-8 °C

After reconstitution: **8 hours at 18 – 25 °C**
1 week at 2 – 8 °C

Warnings and Precautions

The blood donations used for production are tested with CE-marked test kits and found to be non-reactive for HBsAg, anti-HIV-1/2 und anti-HCV. In addition HIV and HCV is tested by PCR. Despite of that the danger of infection can not be excluded with certainty. The control materials have to be treated with the same care as patient samples.³⁾

Waste Management

Dispose of any discarded materials in accordance with the requirements of local regulations authorities.

Limitations

The products should not be used past the expiration date.

The products are not intended for use as a calibrator.

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DiaServe Laboratories GmbH
Scil Diagnostics GmbH