6081 UriTrol, Level 1
6082 UriTrol, Level 2

10 x 5.0 ml
10 x 5.0 ml

For In Vitro Diagnostic Use
In Vitro Diagnostika

Intended Use

UriTrol, Level 1/2 are used as assayed urine controls for Internal Quality Control (IQC) in medical laboratories¹⁾.

Principle

The use of control materials is an objective assessment to monitor precision and accuracy. Two levels of control are available to allow performance monitoring.

Reagent

UriTrol, Level 1/2 are lyophilized control materials based on human urine for the following analytes:

Albumin, Amylase, Calcium, Chloride, Creatinine, Glucose, Phosphorous, Magnesium, Potassium, Sodium, Protein, Urea, Uric Acid.

Assigned Values and Control Ranges

Values and Control Ranges are based on replicate measurements in several laboratories and are specific for the lots of UriTrol, Level 1/2.

Assigned values are only valid for the indicated method and may differ due to various test kits, calibrators and system applications.

Procedure

UriTrol, Level 1/2 should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument or reagent kit being used.

Reconstitution

Remove screw cap from the vial and pipette exactly 5.0 ml dest. water to the lyophilisate.

Keep the vial for about 30 min at rest until the lyophilisate has completely dissolved.

Swing carefully- do not shake! Avoid developing of foam.

Storage and Stability

UriTrol, Level 1/2 are stable until the expiration date when stored unopened at +2 °/+8 °C.

After reconstitution: **8 hours at +18 °/+25 °C**
2 days at +2 °/+8 °C
1 month at - 20 °C

Warnings and Precautions

The blood of the urine donors were tested with CE-marked test kits and found to be non-reactive for HBsAg, anti-HIV-1/2 und anti-HCV. In addition HIV and HCV were tested by PCR. Despite of that the danger of infection from biological material can not be excluded with certainty. The control materials have to be treated with the same care as patient samples²⁾.

Waste Management

Dispose of any discarded materials in accordance with the requirements of local regulations authorities.

Limitations

The products should not be used past the expiration date.

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DiaServe Laboratories GmbH
Scil Diagnostics GmbH